

The Afinion Platform

Product Tailored For Physician Office Testing



Simple operation – no risk of incorrect results (CLIA-waived)



Fast – results during consultation



Reliable – lab-equivalent precision and accuracy

High standards that benefit all markets

Best Performance In POC Testing



Simplest test procedures



Automated sample processing



Most rapid HbA1c result



Best-in-class HbA1c performance




Bi/directional connectivity



Panel of multiple tests during the consultation

HbA1c, ACR, Lipid Panel, CRP

Abbott Afinion System




Afinion 2 key features



- Quality
- Menu
- Ease of use
- Connectivity
- Reliable
- Quiet
- Speed

Afinion Assays



Simple 3-step procedure



1. Collect a sample from the finger, venous tube or urine cup using the integrated sampling device
2. Place the sampling device back into the test cartridge
3. Place the cartridge in the Analyzer and close the lid. Processing will start automatically. The test result will be displayed in minutes.

Afinion 2 – Simply more efficient.

